

Venous Blood Specimen Collection


- ✓ Screw-on needle hub compatible with all available
- ✓ Precision lancet point grinding, silicone treated tip
- ✓ Easy-to-open packaging and manipulation
- ✓ Comfortable collection experience
- ✓ Sterile, individually packed
- ✓ Intended for blood collection with vacuum tubes
- ✓ Made of polypropylene hub and stainless steel tube
- ✓ Thin wall tube material
- ✓ Color coded needle hub

Description	Cat. No.	Color Code	Size G/”	Diameter	Length
Blood specimen Collection Needle	CY1162425		24G x 1”	0.55 mm	26 mm
	CY1162325		23G x 1”	0.60 mm	32 mm
	CY1162225		22G x 1”	0.70 mm	35 mm
	CY1162125		21G x 1”	0.80 mm	35 mm
	CY1162025		20G x 1”	0.90 mm	35 mm
	CY1161825		18G x 1”	1.20 mm	35 mm

*** Packing info**

Product	Needle Size	Primary	Secondary	Packing Size	G.W.	N.W.
Blood specimen Collection Needle	18G	100 pcs / box	50 box / Ctn	45*24*44 cm	12.4 Kg	11.6 Kg
	20G,21,22 & 23G	100 pcs / box	50 box / Ctn	45*24*44 cm	13.2 Kg	12.4 Kg


Blood Collection Needle Holder

- Single Use Device: ✓ Eliminates needle removal and holder reuse
 Total Needle Protection: ✓ Front and back end maximizes protection coverage
 No Technique Change: ✓ Promotes ease of use.
 Safety By Design: ✓ The safe and easy solution to your need for a single use device

Product	Cat. No.	Size	Qty.per bag/case	Packing Size	G.W.	N.W.
Needle holder	CY133A	20*55 mm	250 / 20	65*40*37 cm	15.0 Kg	14.0 Kg