

200 uL Pipette tips

Ensure accurate, precise pipetting.

Manufactured using high quality raw materials and the latest molding techniques. Designed using our knowledge of liquid handling. Our modern automated production facilities, and our expertise, enable us to produce uniformly smooth, hydrophobic surfaces on the inside and outside of our tips, to prevent liquid retention.

- ✓ Available sterile or nonsterile
- ✓ Packaging options include racks, bags, and boxes.
- ✓ Certified to be free from DNA, DNase, RNase and endotoxin
- ✓ Tips come vacuum sealed in irradiation-sterilized hinged tip racks
- ✓ Smooth, uniform surface prevents liquid retention

Order Information

200 uL with filter, DNase/Rnase-free, Non-Pyrogenic, Sterile/Non-Sterile available

| Cat. No. | Product | Filter | Low Retention | Capacity | Color | Packing | Qty. per bag/case |
|-----------|---------------|--------|---------------|----------|-----------|----------------|-------------------|
| PPT151200 | Filtered Tips | ● | | 200 µL | Colorless | Resealable bag | 1000 / 20 |
| PPT153200 | Filtered Tips | ● | | 200 µL | Colorless | Racked | 96 / 50 |
| PMT111200 | Filtered Tips | ● | ● | 200 µL | Colorless | Resealable bag | 1000 / 20 |
| PMT231200 | Filtered Tips | ● | ● | 200 µL | Colorless | Racked | 96 / 50 |





200 uL, DNase/Rnase-free, Non-Pyrogenic, Sterile/Non-Sterile available

| Cat. No. | Product | Filter | Low Retention | Capacity | Color | Packing | Qty. per bag/case |
|-----------|-------------------|--------|---------------|----------|--------|----------------|-------------------|
| PPT001200 | Pipette Tips | | | 200 µL | Yellow | Resealable bag | 1000 / 20 |
| PPT011200 | Pipette Tips | | | 200 µL | Yellow | Racked | 96 / 50 |
| PMT011200 | Low Retention Tip | | ● | 200 µL | Yellow | Resealable bag | 1000 / 20 |
| PMT251200 | Low Retention Tip | | ● | 200 µL | Yellow | Racked | 96 / 50 |



Technical Specification



200 uL filtered Tips