

## 20 uL Pipette tips

Ensure accurate, precise pipetting.

Manufactured using high quality raw materials and the latest molding techniques. Designed using our knowledge of liquid handling. Our modern automated production facilities, and our expertise, enable us to produce uniformly smooth, hydrophobic surfaces on the inside and outside of our tips, to prevent liquid retention.

- ✓ Available sterile or nonsterile
- ✓ Packaging options include racks, bags, and boxes.
- ✓ Certified to be free from DNA, DNase, RNase and endotoxin
- ✓ Tips come vacuum sealed in irradiation-sterilized hinged tip racks
- ✓ Smooth, uniform surface prevents liquid retention


### Order Information

20 uL with filter, DNase/Rnase-free, Non-Pyrogenic, Sterile/Non-Sterile available

Cat. No.	Product	Filter	Low Retention	Capacity	Color	Packing	Qty. per bag/case
PPT101020	Filtered Tips	●		20 µL	Yellow	Resealable bag	1000 / 20
PPT111020	Filtered Tips	●		20 µL	Yellow	Racked	96 / 50
PMT111020	Filtered Tips	●	●	20 µL	Yellow	Resealable bag	1000 / 20
PMT252020	Filtered Tips	●	●	20 µL	Yellow	Racked	96 / 50



### Technical Specification



**20 uL filtered Tips**