

Cell Culture Dishes

Cell and Tissue Culture Dishes are ideal containers for cell growth and yields aim on little and medium input volume, and also useful in sample separation, pre-treatment, storage and so on.

- ✓ Flat bottom uniform wall thickness ensures distortion-free bottom
- ✓ Dish surface is smooth and free from striation to maximize usable area for growth
- ✓ Lids with several little chimps to shape vents are available for very effective gas exchange
- ✓ Package in durable zip resealable self-standing plastic bags that allow dishes to remain upright and reduce contamination
- ✓ The rim on upper side of the lid mates with the dish brim for easy and secure stack
- ✓ *DNase/RNase-free and Non-pyrogenic

Order Information



Cat. No.	Product	Size	Cell Growth Area	Surface Treated	Qty. per bag/case
TCD000035	Cell Culture Dishes	Ø35 × 13 mm	8.8 cm ²		10 / 96
TCD000060	Cell Culture Dishes	Ø60 × 16 mm	21.5 cm ²		10 / 60
TCD000100	Cell Culture Dishes	Ø100 × 17 mm	55.1 cm ²		10 / 30
TCD000150	Cell Culture Dishes	Ø150 × 21 mm	143.0 cm ²		1 / 120
TCD010035	Cell Culture Dishes	Ø35 × 13 mm	8.8 cm ²	●	10 / 96
TCD010060	Cell Culture Dishes	Ø60 × 16 mm	21.5 cm ²	●	10 / 60
TCD010100	Cell Culture Dishes	Ø100 × 17 mm	55.1 cm ²	●	10 / 30
TCD010150	Cell Culture Dishes	Ø150 × 21 mm	143.0 cm ²	●	1 / 120

Bacteriological Petri Dishes

Single-use Petri dishes for microbiological procedures, as well suspension culture or plant cell culture

- ✓ Non-cytotoxic virgin polystyrene, excellent optical clarity
- ✓ Engineered for optimum flatness to provide uniform agar thickness
- ✓ Ø90 and 150 mm Petri dishes are manufactured in accordance with standard ISO 24998:2008, ideal for use with most automatic plate pourers.



Cat. No.	Product	Size	Room	Vents	Packing	Qty. per bag/case
BPD000035	Bacteriological Petri Dishes	Ø35 × 15 mm	1		polybag	10 / 200
BPD000060	Bacteriological Petri Dishes	Ø60 × 15 mm	1	●●●	polybag	10 / 200
BPD000070	Bacteriological Petri Dishes	Ø70 × 15 mm	1		polybag	10 / 100
BPD000390	Bacteriological Petri Dishes	Ø90 × 15 mm	1	●●●	polybag	10 / 50
BPD000090	Bacteriological Petri Dishes	Ø90 × 15 mm	1		polybag	10 / 50
BPD000150	Bacteriological Petri Dishes	Ø150 × 15 mm	1	●●●	polybag	10 / 20
BPD000092	Bacteriological Petri Dishes	Ø90 × 15 mm	2	●●●	polybag	10 / 50
BPD000093	Bacteriological Petri Dishes	Ø90 × 15 mm	3	●●●	polybag	10 / 50
BPD000094	Bacteriological Petri Dishes	Ø90 × 15 mm	4		polybag	10 / 50
BPD000130	Petri Dishes - Square	130×130×15 mm	1		polybag	10 / 30

Inoculation Loops



Cat. No.	Product	Style	Inoculation Aerea	Sterile	Qty. per bag/case
INL000110	Inoculation Loops	Square	1 µL+ 10 µL	●	5 / 1000
INL000201	Inoculation Loops	Needle	1 µL	●	5 / 1000
INL000210	Inoculation Loops	Needle	10 µL	●	5 / 1000
INL000301	Inoculation Loops	Sling	1 µL	●	5 / 1000
INL000310	Inoculation Loops	Sling	10 µL	●	5 / 1000