

PCR 8 Tubes Strip




- ✓ Strip of 8 tubes with Flat caps
- ✓ Angled caps that provide efficient tube strip positioning in racks and thermal cycler blocks
- ✓ Increased labeling area: dual end tabs for proper labeling and easy handling
- ✓ Graduated 1 to 8 number markers on each tube for visual verification
- ✓ Graduated 20 μ L measuring markers on each tube for visual volume verification
- ✓ Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.

Order Information



Cat. No.	Product	Capacity (Metric)	Closure Type	Packing	Qty. per bag/case
PCR100200	PCR 8 Tubes Strip	0.2 mL	Flat	Resealable bag	125 / 10

Technical Specification



PCR 8 Tubes Strip