

Centrifuge tube 15 mL & 50 mL

Made from high-purity resins and molded using state-of-the-art processes; safe alternative to glass without sacrificing accuracy meets FDA 21 CFR 177.1520 & USP Class VI.

- ✓ Flat Cap
- ✓ Non-Cytotoxic
- ✓ RNase/DNase-free
- ✓ Biologically inert: USP biological reactivity class VI
- ✓ High level of sterility: SAL 10⁻⁶
- ✓ Leak-proof to protect your samples
- ✓ Max RCF up to 12,000 x g
- ✓ Graduations every 5 mL increments.
- ✓ Nonpyrogenic: Endotoxin tested (0.5 EU/mL)

Order Information




Cat. No.	Product	Max. RCF	Capacity	Bottom	Sterile	Packing	Qty. per bag/case
CFT000015	Centrifuge tube	12,000	15 mL	Conical	●	Resealable bag	50 / 10
CFT001015	Centrifuge tube	12,000	15 mL	Conical	●	Resealable bag	50 / 10
CFT000020	Centrifuge tube	12,000	50 mL	Conical	●	Resealable bag	25 / 20
CFT001020	Centrifuge tube	12,000	50 mL	Conical	●	Resealable bag	25 / 20

Technical Specification



Centrifuge tube 50 mL, Flat Cap



Centrifuge tube 15 mL, Flat Cap