

Cryogenic-Storage Tube

Small, capped vials designed for storing liquid samples at ultra-low temperatures. Vials come in various volume capacities and designs and are suitable for general cold storage or vapor state liquid nitrogen storage. Sterile and non-sterile formats are available.

- ✓ Large writing areas accommodate more written sample information and graduations for easy visual volume assessment
- ✓ Noncytotoxic raw materials for added sample security and tubes certified free of DNase and RNase for genomic applications
- ✓ Certified SAL of 10⁻⁶ for the highest sterility assurance available
- ✓ Suitable for storage from general cold storage refrigeration (+4°C) down to the vapor phase of liquid nitrogen
- ✓ Self standing tubes for ease of use on the benchtop
- ✓ Storage temperature range -196°C ~ 121°C


Order Information

Cryogenic-Storage Tube with externally threaded cap


Cat. No.	Product	Capacity	Thread Style	Sterile	Packing	Qty. per bag/case
CST511105	Cryogenic-Storage Tube	0.5 mL	Externally	●	Resealable bag	50 / 20
CST511205	Cryogenic-Storage Tube	0.5 mL	Externally		Resealable bag	50 / 20
CST511110	Cryogenic-Storage Tube	1.0 mL	Externally	●	Resealable bag	50 / 20
CST511210	Cryogenic-Storage Tube	1.0 mL	Externally		Resealable bag	50 / 20
CST511118	Cryogenic-Storage Tube	1.8 mL	Externally	●	Resealable bag	50 / 20
CST511218	Cryogenic-Storage Tube	1.8 mL	Externally		Resealable bag	50 / 20
CST511136	Cryogenic-Storage Tube	3.6 mL	Externally	●	Resealable bag	50 / 20
CST511236	Cryogenic-Storage Tube	3.6 mL	Externally		Resealable bag	50 / 20
CST511145	Cryogenic-Storage Tube	4.5 mL	Externally	●	Resealable bag	50 / 20
CST511245	Cryogenic-Storage Tube	4.5 mL	Externally		Resealable bag	50 / 20
CST512110	Cryogenic-Storage Tube	0.5 mL	Internally	●	Resealable bag	50 / 20
CST512210	Cryogenic-Storage Tube	0.5 mL	Internally		Resealable bag	50 / 20
CST512118	Cryogenic-Storage Tube	1.0 mL	Internally	●	Resealable bag	50 / 20
CST512218	Cryogenic-Storage Tube	1.0 mL	Internally		Resealable bag	50 / 20
CST512136	Cryogenic-Storage Tube	1.8 mL	Internally	●	Resealable bag	50 / 20
CST512236	Cryogenic-Storage Tube	1.8 mL	Internally		Resealable bag	50 / 20
CST512145	Cryogenic-Storage Tube	3.6 mL	Internally	●	Resealable bag	50 / 20
CST512245	Cryogenic-Storage Tube	3.6 mL	Internally		Resealable bag	50 / 20



Technical Specification



1.8 mL Externally-Thread Tube



1.8 mL Internally-Thread Tube